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J1 Pin	Function	A/D Pin Connection
1	Analog Channel 0	37
2	Analog Channel 1	36
3	Analog Channel 2	35
4	Analog Channel 3	34
5	Analog Ground	19
6	Reserved	
7	N/C	
8	N/C	
9	Digital Ground	7
10	A/D Trigger	25
11	Master Clock	20
12	Digital Ground	Not in Cable
13	+5VDC	1
14	Analog Channel 4	33
15	Analog Channel 5	32
16	Analog Channel 6	31
17	Analog Channel 7	30
18	Analog Ground	Not in Cable
19	N/C	
20	Continuous Strobe In	8 (or use internal jumpers)
21	External Software Trigger Out (DO3)	5
22	Spectrometer Mode Input S1	4
23	Integration Time Clock In	2
24	Strobe Enable, Spectrometer Mode Input S0	23
25	Enable Read In	3

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Icon	Function	Description
	Open	Opens a new spectral window.
	Open a Processed Spectrum	Opens a processed spectrum and displays the data in the spectral window.
	Save Processed Spectrum	Saves the processed spectra to disk.
	Copy	Copies the current spectra to the clipboard.
	Print	Prints the currently displayed spectra.
	Help	Opens the OOIBase32 help system.
	Store Global Dark	Take a dark spectrum used by all spectral windows.
	Store Global Reference	Take a reference spectrum used by all spectral windows.
	Global Snapshot	Take a snapshot and freeze data acquisition on all spectral windows.
	Global Emergency Reset	Reset the acquisition parameters for all spectral windows.
	Kick Start	Restart the acquisition loop without resetting any acquisition parameters.
	Store Dark	This option stores a dark spectrum for all enabled spectrometer channels in a spectral window. Block the light path to the sample, and then take the dark spectrum.
	Store Reference	This option stores a reference spectrum for all enabled spectrometer channels in a spectral window. Take a reference spectrum with the light source on and a blank in the sampling region.
	Snapshot	This option halts data acquisition and takes a snapshot of the activity in the spectral window.
	Single Exposure	This option reactivates data acquisition, and acquires and displays a single scan. It is only active when OOIBase32 is in Snapshot mode.
	Configure Data Acquisition	This option opens the Configure Data Acquisition dialog box. This dialog box allows you to configure aspects of the data acquisition process.
	Emergency Reset	This option resets all acquisition parameters for the active spectral window.
	Toggle Cursor	Enables or disables the display of a vertical cursor for the spectral window.
	Cursor Peak Left	Moves the cursor to the next left peak.
	Cursor Big Left	Moves the cursor 25 pixels to the left.
	Cursor Left	Moves the cursor 1 pixel to the left.
	Cursor Right	Moves the cursor 1 pixel to the right.
	Cursor Big Right	Moves the cursor 25 pixels to the right.
	Cursor Peak Right	Moves the cursor to the next right peak.

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Icon	Function	Description
	Configure Cursor	Opens the Configure Cursor dialog box.
	Autoscale	Autoscale the graph to fit the spectral window.
	Set Scale	Set the scale of the graph.
	Unscale	Restores the graph to the default scale.
	Subtract Dark Spectrum	This command switches the current spectral window into Scope Mode, and subtracts the stored dark spectra from each spectrometer channel before OOIBase32 displays it.
	Scope Mode	This command switches the current spectral window into Scope Mode.
	Absorbance Mode	This command switches the current window into Absorbance Mode.
	Transmission Mode	This command switches the current window into Transmission Mode.
	Relative Irradiance Mode	This command switches the current window into Relative Irradiance Mode.
	Specular Reflection Mode	This command switches the current window into Specular Reflection Mode.
	Script-defined Custom Mode	This mode is only available in OOIBase32 Platinum version.
	Configure Spectrometer	Opens the Spectrometer Configuration dialog box.
	Activate Time Acquisition Mode	Places OOIBase32 in Time Acquisition Mode.
	Start Time Acquisition	Starts the time acquisition process.
	Pause Time Acquisition	Pauses the time acquisition process.
	Stop Time Acquisition	Stops the time acquisition process.
	Suspend Graph Updates	Suspends the graph display during a time acquisition process.

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Spectrometer Configuration

Wavelength Calibration | A/D Interface | Reference Monitoring | Stray Light Correction | Detector Linearity

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Option	Description
Spectrometer Channel	Specifies the spectrometer channel for which these modifications will apply
Spectrometer Serial Number	Serial number of the spectrometer on the channel selected in the Spectrometer Channel section
Channel Enabled	Specifies whether or not the channel specified in the Spectrometer Channel section will acquire data
First Coefficient	First wavelength coefficient (provided on the Wavelength Calibration Data Sheet that came with the spectrometer). Applies to the spectrometer selected in the Spectrometer Channel section
Second Coefficient	Second wavelength coefficient (provided on the Wavelength Calibration Data Sheet that came with the spectrometer). Applies to the spectrometer selected in the Spectrometer Channel section
Third Coefficient	Third wavelength coefficient (provided on the Wavelength Calibration Data Sheet that came with the spectrometer). Applies to the spectrometer selected in the Spectrometer Channel section. If your Wavelength Calibration Data Sheet does not contain a third coefficient, enter 0 here
Intercept	Specifies the wavelength intercept of the spectrometer selected in the Spectrometer Channel section

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Copyright © by Chiang Mai University
All rights reserved