

LIST OF TABLES

Table		Page
1	Important parameters in various steps of sol – gel process.....	10
2	Classification of infrared radiation.....	24
3	Unit cell parameters of crystal structure.....	38
4	XPS binding energy (BE) values of references.....	45
5	Sample codes of the catalysts prepared by the Pechini method according to heating conditions.....	58
6	Sample codes of the catalysts prepared by Schiff base complex method according to heating conditions.....	60
7	Sample codes of the catalysts prepared by coprecipitation according to heating conditions.....	65
8	Instrument broadening (FWHM, $\beta_{\text{reference}}$).....	68
9	The energy scanning ranges (eV), energy step and number of scan of the investigated elements.....	70
10	Reaction conditions for the catalytic oxidation of toluene on LaCoO ₃ catalysts.....	76
11	FTIR spectra assignment of samples prepared by the Pechini method.....	80
12	Diffraction data of LaCoO ₃ reference and LaCoO ₃ prepared by the Pechini method.....	85
13	Calculated values of crystallite sizes of samples prepared by the Pechini method.....	88
14	FTIR spectra assignment of samples prepared by Schiff base complex method.....	91
15	Diffraction data of LaCoO ₃ reference and LaCoO ₃ prepared by Schiff base complex method.....	96
16	Calculated values of crystallite sizes of samples prepared by Schiff base complex.....	99

LIST OF TABLES (continued)

Table		Page
17	FTIR spectra assignment of samples prepared by coprecipitation.....	100
18	Diffraction data of LaCoO ₃ reference and LaCoO ₃ prepared by coprecipitation.....	103
19	Diffraction data of La(OH) ₃ reference and La(OH) ₃ impurity in LaCoO ₃ prepared by coprecipitation.....	104
20	Calculated values of crystallite sizes of samples prepared by coprecipitation.....	105
21	Binding energies of La 3d of the prepared LaCoO ₃	109
22	Binding energies of O 1s of the prepared LaCoO ₃	112
23	Binding energies of Co 2p of the prepared LaCoO ₃	116
24	The XRF results of the prepared LaCoO ₃	117
25	Surface area of the prepared LaCoO ₃	118
26	DRIFTS spectra assignment of LC-PN773-08 catalyst.....	120
27	The time used to convert toluene to 90% related to that of zero time.....	124
Appendix Table		
A1	FTIR spectra assignment of Co(NO) ₃ · 6H ₂ O calcined at 773 K....	140
A2	The JCPDS file of LaCoO ₃ with the number: 25-1060.....	141
B1	Binding energies of C 1s for LaCoO ₃ prepared by the Pechini method.....	144
B2	Binding energies of C 1s for LaCoO ₃ prepared by the Schiff base complex method.....	146
C1	The flow rate of toluene at various pump drive speeds.....	147
D1	Sample codes according to the preparation temperature.....	152

LIST OF TABLES (continued)

Appendix Table		Page
D2	The effect of the calcination temperature on the catalytic activity	156